


www.medaid.co.uk



The MAI operating table



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Tel: 01234 930 394
tim@medaid.co.uk
www.medaid.co.uk

 /medicalaidinternationalcic



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments



INTRODUCING The MAI operating table

A purpose built operating table designed for the challenges of low resource environments

- > Easy to transport, can be taken apart for outreach or disaster relief work, simply remove four bolts
- > Supplied with essential attachments (arm boards X 2, arm table, lithotomy poles)
- > Wide range of additional attachments available
- > Wide variety of configurations, suitable for all types of surgery
- > CE marked, high quality European manufacture
- > Manual operation, no electricity required
- > Stainless steel
- > Value for money
- > Variety of bases available



Simple to take apart for transport



Arm table for hand and arm operations



Arm board for easy IV access



A variety of lithotomy poles are available for gynaecological/obstetric and urological procedures



Trendelenburg to 45°, ideal for Fistula procedures



Shoulder supports for use when steep Trendelenburg is required



Sitting position for shoulder procedures



Fracture table for lower limb trauma



Option to remove or angle individual legs at bottom of table

