

Paediatric Thoracotomy Set List

Quantity	Item
1	Weinstein Instrument Holder
2	Halstead Mosquito Artery Forceps Box Joint 130mm Cvd
4	Crile Artery Forceps Box Joint 140mm CVD
2	Crile Artery Forceps Box Joint 140mm STR
2	Rochester Pean Artery Forceps Box Joint 130mm CVD
2	Rochester Pean Artery Forceps Box Joint 140mm CVD
2	Rochester Pean Artery Forceps Box Joint 160mm CVD
2	Allis Tissue Forceps 4 Into 5 Teeth Box Joint 150mm
2	Allis Tissue Forceps 5 Into 6 Teeth Box Joint 150mm
2	Allis Tissue Forceps 5 Into 6 Teeth Box Joint 180mm
2	Babcock Tissue Forceps Box Joint 160mm
2	Babcock Tissue Forceps Box Joint 180mm
2	Rochester Oschner Artery Forceps Box Joint 1 Into 2 Teeth 160mm STR
2	Schmidt Artery Forceps Box Joint 190mm CVD
1	Bridge Forceps Curved
2	Baby Mixer Artery Forceps Box Joint 140mm CVD
2	Mixer O'Shaughnessy Artery Forceps Box Joint 190mm CVD
1	Mason Needle Holder TC Jaws and Box Joint 265mm (Masson) STR
1	Mayo Hegar Needle Holder TC Jaws Box Joint 140mm
1	Mayo Hegar Needle Holder TC Jaws Box Joint 150mm
1	Mayo Hegar Needle Holder TC Jaws Box Joint 180mm
1	Mayo Scissors TC Jaws 140mm STR
1	Metzenbaum Scissors TC Jaws 165mm CVD
1	Metzenbaum Scissors TC Jaws 180mm CVD
1	Metzenbaum Scissors TC Jaws 200mm CVD
1	Mayo Scissors TC Jaws 180mm STR
1	Debaquey Tangent Vascular Clamp
1	Debaquey Tangent Vascular Clamp
1	Debaquey Tangent Vascular Clamp
2	Duval Lung Forceps
1	Foerster Dressing + Sponge Holding Forceps Serrated Jaws + Box Joint 180mm STR
2	Lorna (Edna) Towel Clip with Box Joint 100mm
2	Backhaus Towel Clip Box Joint 115mm CVD

Quantity	Item
2	Scalpel Handle No.3
2	Continental Pattern Dissecting Forceps Serrated Jaws and Narrow End 130mm
2	Continental Pattern Dissecting Forceps Serrated Jaws and Narrow End 140mm
2	Continental Pattern Dissecting Forceps Serrated Jaws and Narrow End 180mm
2	Lane Tissue Forceps 1 Into 2 Teeth Box Joint 130mm
2	Lane Tissue Forceps 1 Into 2 Teeth Box Joint 195mm
2	Debakey Tissue Dissecting Forceps Atraumatic Jaws 150mm
2	Debakey Tissue Dissecting Forceps Atraumatic Jaws 200mm
2	Potts Smith Dissecting Forceps TC Jaws Narrow 180mm STR
2	Russian Pattern Dissecting Forceps Oval Toothed End 150mm STR
2	Adson Dissecting Forceps 1 Into 2 Teeth 130mm
1	US Army Navy Retractor
1	Cushing Vein Retractor
2	Richardson Retractor 28mm Wide X 20mm Deep 240mm
2	Richardson Retractor 36mm Wide X 28mm Deep 240mm
2	Kelly Retractor Small Size 38mm End Of Blade 190mm
1	Balfour Baby Abdominal Retractor with Centre Blade and Fenestrated End Blades
1	Finochietto Rib Spreader Small Size with 30mm x 45mm Blades and 165mm Spread
1	Finochietto Rib Spreader Baby Size with 18mm x 21mm Blades and 75mm Spread
1	Finochietto Rib Spreader Medium Size, 45mm x 65mm Blades and 200mm Spread
1	Bailey Rib Spreader 160mm
1	Bailey Rib Contractor (Short Claws)
1	Bailey Rib Contractor (Long Claws)
1	Allison Retractor Infant Size 280mm
1	Allison Retractor Child Size 280mm